



Recombinant Acanthamoeba polyphaga mimivirus Putative SWIB domain-containing protein R508 (MIMI_R508)

Product Code	CSB-BP716480ADAZ
Abbreviation	MIMI_R508
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UQ74
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MSKRVTSSKK SKIMREDSSEN DSAKNLSKTS KSVSKSKEAG SKRTKQNSSK SSKDLKTTKN IGGSKSSRTY NSEGSKNIKK TSSSKDSKVI KKNKQKVESS DSEKHSENKS HKKSSKSSSI SRKKPIKKDT RKINIELVSD SDNNSDSELN GVDTTGDDII YLKGEDNNDK QQNKKSPKKL LNEKKISSES FDDKLNELRE ELRENYIRQK KIMNDIKDLA NLHKKDLRLA VKSGNRLNSG KHSGFNKPQT VPQSLKDLLK INDDVLSRSK VTELMYKYFT DNKMYNSKTK REIIPNSKIK KLFGMKEGDI ITFYNMQTWL KKVYNENQNN EKLEE
Source	Baculovirus
Target Names	MIMI_R508
Protein Names	Recommended name: Putative SWIB domain-containing protein R508
Expression Region	1-335
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.