



Recombinant Acanthamoeba polyphaga mimivirus

Putative ankyrin repeat protein R911 (MIMI_R911)

Product Code	CSB-EP716515ADAZ-B
Abbreviation	MIMI_R911
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UR04
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MYDTLPPELW VQIVDYSGEI SLLLTNTSFF ELFNLINIKT DII EYVMEND LIDILKYL SL LKGLGHP IIV DKNILNIKT L NKYLIKNCGE NQLEIVKFLV SLGADIRAGN DYAVGLSSQN GHLEVVKYLV NQGSDIRAEN DYAVRWASGN GHLEVVKYLV SQGANIRADN DHAIGLASYY GYLEVVKYLV SQGADIRSDN DYAVRMASRN GHIEVVEYLV SQGANIRSDN DYAVRLASQN GHLEVVKYLV SQGADIKSDN DYAVRLASQN GHLEVVEYLV TQGTNIRVNN NYAVEWASKN GNLEVVKYLI SQGADIADN NFAVRWASRN GHLEVVKYLV SLGADIKSDN DYAVRWASGN GHLEVVKYLV SQGSDIRVEN DYAVRWASRN GHFDVIKYL V SQGADIRSDN DYAVKWASEN GHLEVVKFLV SLGADIKAED DYAVRWASEK GHLEVVEYLV SQGAVLS
Source	E.coli
Target Names	MIMI_R911
Protein Names	Recommended name: Putative ankyrin repeat protein R911
Expression Region	1-467
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.