



# Recombinant Acanthamoeba polyphaga mimivirus Uncharacterized N-acetyltransferase R584 (MIMI\_R584)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP716521ADAZ-B                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Abbreviation</b>      | MIMI_R584                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                       |
| <b>Uniprot No.</b>       | Q5UR52                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Immunogen Species</b> | Acanthamoeba polyphaga mimivirus (APMV)                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Sequence</b>          | MPKSSSNNPC PPGEILREGY YRHGYERRGF KRSSGTYIPP TEVSEAYVPP<br>TCIIDRGKPG KGPRILPKPG NDLHLSWYGY AVHEPERDRH RALVEASQEN<br>DPLEVLRRRLN LLRNYQAYPD VKDIMSQDVK FMSDFYAKYK ERNQEANRFG<br>RSNSRSRSNS RSKSSRSRSN NRSKSSRSSS TQSKSNNRSN SRSNSKTSRR<br>LKKLPENDLD LDQDGGNRYG GDPLTSITSD TTGNQILETN IKLSKETECs<br>EGKCQVFNKV YESHTINGKK IVFETLDQND SENVLSLDKL YLDSQQTIDR<br>VRQNLSTNKG YIIGIKSDNK LEGYCWYKEI SNNEVKIFWF CANKGYGTAL<br>YTFMEKYFKL NNYSRIVIDV SMEGSYAVRR INFWNHQGFK TYQVKTDNHK<br>IHMEKDI |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Target Names</b>      | MIMI_R584                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Protein Names</b>     | Recommended name: Uncharacterized N-acetyltransferase R584 EC= 2.3.1.-                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Expression Region</b> | 1-407                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                              |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.                                                                                                                                                                                                                                                                                                                                               |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.