



Recombinant Acanthamoeba polyphaga mimivirus Uncharacterized N-acetyltransferase R584 (MIMI_R584)

Product Code	CSB-MP716521ADAZ
Abbreviation	MIMI_R584
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UR52
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MPKSSSNNPC PPGEILREGY YRHGYERRGF KRSSGTYIPP TEVSEAYVPP TCIIDRGKPG KGPRILPKPG NDLHLSWYGY AVHEPERDRH RALVEASQEN DPLEVLRRRLN LLRNYQAYPD VKDIMSQDVK FMSDFYAKYK ERNQEANRFG RSNSRSRSNS RSKSSRSRSN NRSKSSRSSS TQSKSNNRSN SRSNSKTSRR LKKLPENDLD LDQDGGNRYG GDPLTSITSD TTGNQILETN IKLSKETECs EGKCQVFNKV YESHTINGKK IVFETLDQND SENVLSLDKL YLDSQQTIDR VRQNLSTNKG YIIGIKSDNK LEGYCWYKEI SNNEVKIFWF CANKGYGTAL YTFMEKYFKL NNYSRIVIDV SMEGSYAVRR INFWNHQGFK TYQVKTDNHK IHMEKDI
Source	Mammalian cell
Target Names	MIMI_R584
Protein Names	Recommended name: Uncharacterized N-acetyltransferase R584 EC= 2.3.1.-
Expression Region	1-407
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.