



Recombinant Haloarcula marismortui Probable DNA primase small subunit (priA)

Product Code	CSB-MP716555HTJ
Abbreviation	priA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5V554
Product Type	Recombinant Protein
Immunogen Species	Haloarcula marismortui (strain ATCC 43049 / DSM 3752 / JCM 8966 / VKM B-1809) (Halobacterium marismortui)
Purity	≥85% (SDS-PAGE)
Sequence	MEERTRAYLR GRFGDHYRQA SVTPPPAANE REWGFIPWTE GPGETMVRHR SLLDLGEIED FLGRRKPRHV YFSAGRYDEP SASTMSDKGW RSSDLVFDLD ADHLPSVVLG EDSYAEMLA CKDALRRLLD FLEDDFGFDD LTIVFSGGRG YHVHVRDERI RHLERDARRE VVDYVRGIGL EFDELVD EES VAGTAGRSSP AQKRTLSTEG GWSARAHHRM LAVVDDLLAM EEADALEQLQ EYDGIGEGKA TAALNAARSN YEQLEAGNID VHPAFYQLAK ILLHEVVAAD NAPIDEPVTT DTNRLIRLPG SLHGGSGLEV QRIDRDDLDA FDPLVDPVPE TFRGHDITVE VTDGGLVELD GDSFTLEAGN QTVPEHVGVF LMARGRAEKG KE
Source	Mammalian cell
Target Names	priS
Protein Names	Recommended name: Probable DNA primase small subunit EC= 2.7.7.-
Expression Region	1-392
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.