



Recombinant Human Protein FAM194B

(FAM194B)

Product Code	CSB-BP716587HU
Abbreviation	FAM194B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5W0A0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSAENNQLSG ASPPHPPTTP QYSTQNLPSE KEDTEVELDE ESLQDESPFS PEGESLEDKE YLEEEEDLEE EEYLGKEEYL KEEEYLGKEE HLEEEYLEK AGYLEEEYI EEEYLGKEG YLEEEYLGK EEHLEEEYLGK GKEGYLEKED YIEEVDYLGK KAYLEEEYLGK GKKSYLEEEK ALEKEENLEE EEALEKEENL DGKENLYKKY LKEPKASYSS QTMLLRDARS PDAGPSQVTT FLTVP LTFAT PSPVSESATE SSELLLTLYR RSQASQTDWC YDRTAVKSLK SKSETEQETT TKLAPEEHVN TKVQQKKEEN VLEFASKENF WDGITDESID KLEVEDLDEN FLNSSYQTVF KTIKEMAAH NELEEDFDIP LTKLLESENK WKLVI MLKKK YEKFKETILR IKRRREAQKL TEMTSFTFHL MSKPTPEKPE TEEIQKPQRV VHHRKKLERD KEWIKKKT VV HQGDGKLILY PNKNVYQILF PDGTGQIHYP SGNLAMLILY AKMKKFTYII LEDSLEGRIR ALINNSGNAT FYDENS DIWL NLSSNLGYF PKDKRQKAWN WWNLNIHVHA PPVQPISLKI NEYIQVQIRS QDKIIFCFTY EQKQICNLG TRYKFVIEV LSEMKKKTIL EAEPGPTAQK IRVLLGKMNR LLNYATTPDL ENFIEAVSIS LMDNKYLKKM LSKLWF
Source	Baculovirus
Target Names	ERICH6B
Protein Names	Recommended name: Protein FAM194B
Expression Region	1-696
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.