



# Recombinant Human Protein FAM194B (FAM194B)

<b>Product Code</b>	CSB-YP716587HU
<b>Abbreviation</b>	FAM194B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5W0A0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSAENNQLSG ASPPHPPTTP QYSTQNLPS E KEDTEVELDE ESLQDESPFS PEGESLEDKE YLEEEEDLEE EEYLGKEEYL KEEEYLGKEE HLEEEYLEK AGYLEEEYI EEEYLGKEG YLEEEYLGK EEHLEEEY L GKEGYLEKED YIEEVDYLGK KAYLEEEY L GKKSYLEEEK ALEKEENLEE EEALEKEENL DGKENLYKKY LKEPKASYSS QTMLLRDARS PDAGPSQVTT FLT VPLTFAT PSPVSESATE SSELLLTLYR RSQASQTDWC YDR TAVKSLK SKSETEQETT TKLAPEEHVN TKVQQKKEEN VLEFASKENF WDGITDESID KLEVEDLDEN FLNSSYQTVF KTIKEMAAH NELEEDFDIP LTKLLESEN R WKLVIMLKK N YEKFKETILR IKRRREAQKL TEMTSFTFHL MSKPTPEKPE TEEIQKPQRV VHHRKKLERD KEWIQKKT VV HQGDGKLILY PNKNVYQILF PDGTGQIHYP SGNLAMLILY AKMKKFTYII LEDSLEGRIR ALINNSGNAT FYDENSDIWL NLSSNLGY YF PKDKRQKAWN WWNLNIHVHA PPVQPISLKI NEYIQVQIRS QDKIIFCFTY EQKQICLN LG TRYK FVIPEV LSEMKKKTIL EAEPGPTAQK IRVLLGKMNR LLNYATTPDL ENFIEAVSIS LMDNKYLKKM LSKLWF
<b>Source</b>	Yeast
<b>Target Names</b>	ERICH6B
<b>Protein Names</b>	Recommended name: Protein FAM194B
<b>Expression Region</b>	1-696
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.