



Recombinant Bothrops jararaca Snake venom serine protease HS114

Product Code	CSB-BP716593BUW
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5W959
Product Type	Recombinant Protein
Immunogen Species	Bothrops jararaca (Jararaca) (Bothrops jajaraca)
Purity	≥85% (SDS-PAGE)
Sequence	VVGGDE CNINEHRSLV AIFNSTGFFC SGILLNQEWV LTASHCDSTN FQMKIGVHSK KTLNQDEQTR NPKEKIFCPN KKNDDALDKD LMLVRLDSPV SDSEHIAPLS LPSSPPSVGS VCRIMGWGSI TPIQKTNPDV PHCANINLLD DAVCRAAYPE LPAEYRTLCA GVPEGGIDTC NGDSGGPLIC NGQFQGIVFY GAHPCGQAPK PGLYTKVIDY NTWIESVIAG NTAATCPP
Source	Baculovirus
Protein Names	Recommended name: Snake venom serine protease HS114 Short name= SVSP EC= 3.4.21.-
Expression Region	25-258
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.