



Recombinant *Bacillus clausii* Ribosomal RNA small subunit methyltransferase H (rsmH)

| | |
|--------------------------|---|
| Product Code | CSB-YP716618BAAC |
| Abbreviation | rsmH |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q5WFG1 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Bacillus clausii</i> (strain KSM-K16) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MFSHITVLKE ESVNGLAVKP GGVYVDCTLG GGGHSERILT ALAGEGHLIA FDQDEAALSF AAEKLSRFKE NITFIRSNFR HIKEELRMRG VDKVDGILFD LGVSSPQLDE ASRGFSYHQD APLDMRMDQS ASLTAREVVN TWPFALKHSI ISRYGEEKFS KQIARKIEAY RLKQPIETT G ELVDIIKEAI PAPARRAGGH PAKRTFQAIR IAVNDELGAF EDALTDGFEL LNEGGRMAVI TFHSLEDRLC KQMFKEKTKL PDLPKGLPII PADQKAPFAL ITKKPIVAND EEIKANNRAR SAKLRIIEKE KQTD |
| Source | Yeast |
| Target Names | rsmH |
| Protein Names | Recommended name: Ribosomal RNA small subunit methyltransferase H EC=2.1.1.199 Alternative name(s): 16S rRNA m(4)C1402 methyltransferase rRNA (cytosine-N(4)-)-methyltransferase RsmH |
| Expression Region | 1-314 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |