



Recombinant Serine hydroxymethyltransferase (glyA)

Product Code	CSB-MP716777SMZ
Abbreviation	glyA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XC65
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
Purity	≥85% (SDS-PAGE)
Sequence	MIFDKGNVED FDKELWDAIH AEEERQEHHI ELIASENMVS KAVMAAQGSV LTNKYAEGYP GNRYYGTEC VDIVETLAIE RAKKLFGAAF ANVQAHSGSQ ANAAAYMALI EAGDTVLGMD LAAGGHLTHG SPVNFSGKTY HFVGYSVDAD TEMLNYEAIL EQAKAVQPKL IVAGASAYS SR SIDFEKFRAI ADHVGAYLMV DMAHIAGLVA AGVHPSPVPY AHIVTSTTHK TLRGPRGGLI LTNDEALAKK INSAVFPGLQ GGP LEHVIAA KAVAFKEALD PAFKDYAQAI IDNTAAMA AV FAQDDRFR LI SGGTDNHVFL VDVTKVIANG KLAQILLDEV NITLNKNAIP FETLSPFKTS GIRIGCAAIT SRGMGVKESQ TIAHLLIKAL VNHNQTVILE EVRQEVRLT DAFPLYKK
Source	Mammalian cell
Target Names	glyA
Protein Names	Recommended name: Serine hydroxymethyltransferase Short name= SHMT Short name= Serine methylase EC= 2.1.2.1
Expression Region	1-418
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.