



Recombinant Phosphoglucosamine mutase (glmM)

Product Code	CSB-EP716780SMZ-B
Abbreviation	glmM
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XCE0
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
Purity	≥85% (SDS-PAGE)
Sequence	MGKYFGTDGV RGEANVELTP ELAFKLGRFG GYVLSQHETE RPKVFFVARDT RISGEMLESA LIAGLLSVGI EVYKLGVLAT PGVSYLVRTE KASAGVMISA SHNPALDNGI KFFGNDGFKL ADDQELEIEA LLDAPEDTLP RPSAEGLGTL VDYPEGLRKY EKFLVTTGTD LSGMTVALDT ANGAASVSAR DVFLDLNAEI AVIGEKPNGI NINDGVGSTH PEQLQELVKE TGADLGLAFD GSDRLIAVD ETGEIVDGD R IMFIIGKYLK EKGLLAHNTI VTTVMSNLGF HKALDKQGIN KAITAVGDRY VVEEMRSSGY NLGGEQSGHV IIMDYNTTGD GQLTAIQLAK VMKETGKSLS ELAAEVTIYP QKLVNIRVEN SMKDRAMEVP AIANIIAKME DEMAGNGRIL VRPSGTEPLL RVMAEAPTA EVDYYVDTIA DVVRTEIGCD N
Source	E.coli
Target Names	glmM
Protein Names	Recommended name: Phosphoglucosamine mutase EC= 5.4.2.10
Expression Region	1-451
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.