



Recombinant Zinc-binding protein AdcA (adcA)

Product Code	CSB-BP716788SMZ
Abbreviation	adcA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XD15
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
Purity	≥85% (SDS-PAGE)
Sequence	EG KVKVVTTFYP VYEFTKGVIG NDGDVSMLMK AGTEPHDFEP STKDIKKIQD ADAFVYMDDN METWVSDVKK SLTSSKVTIV KGTGNMMLLVA GAGHDHHHED ADKKHEHNKH SEEGHNHAFD PHVWLSPYRS ITVVENIRDS LSKAYPEKAE NFKANAATYI EKLKELDKDY TAALSDAKQK SFVTQHAAFG YMALDYGLNQ ISINGVTPDA EPSAKRIATL SKYVKKYGIK YIYFEENASS KVAKTLAKEA GVKAAVLSPL EGLTEKEMKA GQDYFTVMRK NLETLRLTDD VAGKEILPEK DTTKTVYNGY FKDKVVKDRQ LSDWSGSWQS VYPYLQDGTL DQVWDYKAKK SKGKMTAAEY KDYYTTGYKT DVEQIKINGK KKTMTFVRNG EKKTFTYTYA GKEILTPKG NRGVRFMFEA KEPNAGEFKY VQFSDHAIAP EKAHFHLYW GGDSQEKLHK ELEHWPTYYG SDLSGREIAQ EINAH
Source	Baculovirus
Target Names	adcA
Protein Names	Recommended name: Zinc-binding protein AdcA
Expression Region	29-515
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.