



# Recombinant Zinc-binding protein AdcA (adcA)

<b>Product Code</b>	CSB-MP716788SMZ
<b>Abbreviation</b>	adcA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5XD15
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	EG KVKVVTTFYP VYEFTKGVIG NDGDVSMLMK AGTEPHDFEP STKDIKKIQD ADAFVYMDDN METWVSDVKK SLTSSKVTIV KGTGNMMLLVA GAGHDHHHED ADKKHEHNKH SEEGHNHAFD PHVWLSPYRS ITVVENIRDS LSKAYPEKAE NFKANAATYI EKLKELDKDY TAALSDAKQK SFVTQHAAFG YMALDYGLNQ ISINGVTPDA EPSAKRIATL SKYVKKYGIK YIYFEENASS KVAKTLAKEA GVKAAVLSPL EGLTEKEMKA GQDYFTVMRK NLETLRLTTD VAGKEILPEK DTTKTVYNGY FKDKVVKDRQ LSDWSGSWQS VYPYLQDGTL DQVWDYKAKK SKGKMTAAEY KDYYTTGYKT DVEQIKINGK KKTMTFVRNG EKKTFTYTYA GKEILTPKG NRGVRFMFEA KEPNAGEFKY VQFSDHAIAP EKAHFHLYW GGDSQEKLHK ELEHWPTYYG SDLSGREIAQ EINAH
<b>Source</b>	Mammalian cell
<b>Target Names</b>	adcA
<b>Protein Names</b>	Recommended name: Zinc-binding protein AdcA
<b>Expression Region</b>	29-515
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.