



Recombinant Thymidylate kinase (tmk)

Product Code	CSB-BP716795SMZ
Abbreviation	tmk
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q5XDM2
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
Purity	≥85% (SDS-PAGE)
Sequence	MITGKLITVE GPDGAGKTTV LEQLIPLLKQ KVAQEILTTR EPGGVAISEY IRELILDINH TAMDPKTELL LYIAARRQHL VEKVLPALEA GQLVFIDRFI DSSVAYQGAG RGLIKADIQW LNEFATDGLE PDLTLYFDVP SEIGLARINA NQQREVNRLD LETIEIHQRV RKGYLALAKE HPKRIVTIDA TKPLKEVSVV ALEHVLALLL A
Source	Baculovirus
Target Names	tmk
Protein Names	Recommended name: Thymidylate kinase EC= 2.7.4.9 Alternative name(s): dTMP kinase
Expression Region	1-211
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>