



Recombinant Rat Uncharacterized protein C10orf88 homolog

Product Code	CSB-MP716829RA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XI46
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	METAIEDAGL DRGPTLTSSW DAACGALTQS LFLTRTGPR A QDLDFEQLLA PPAQGPDLVS LKSSLSPRDE NPCFIYKCG PNGGEEILSV GILSSARNME VYLGEEYCGT SRGKNVCTVL DNSEHEKILL YKKYLKLESP THACKIKLLS FGEEQCVLIS KVVVHLRPRS AKPSPSSPAL GSRIDLDNIQ TIMESMGSKL SPGAQQLMDM IRFQQQNCLP IRDQLQSVLG TAGHKHLMAL QSSPSSGVVD KASSTPFPR TGLTPSAITE NLKALIDKST QPSGEGNTTN HNECHLMPQN HSLESDLKNA VSSFLPKKAS GSSSVPNSEL LPFLQNLCSQ VNHLRVGHNA RWQENVSKPR EGTVGVPLEE QPVCSYLEKI LSKNMEVMEK KLMKHIDERI YQLQEHIDAK MALLVDLLRS PNSPSPGMPL RHYDSRERLS NGER
Source	Mammalian cell
Protein Names	Recommended name: Uncharacterized protein C10orf88 homolog
Expression Region	1-444
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.