



Recombinant Rat Ubiquilin-like protein (Ubqlnl)

Product Code	CSB-MP716846RA
Abbreviation	Ubqlnl
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XIP4
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MPHIISRSPR MPQSRHPAGV PADGNIPPGV IRVIVKTPGN QNVFTVAEDT SVRQFKEQLS SHFKCQMEQL VLVSMGRLLK DHDTLSQIRGI TDGHTIHVVI KSKHGSRSLT HSFRNLLTNN PCHQDRNPKG NSSTVCQSAG MSETKVESSI LMEPDASKVC TQDPEVGSPE QIEQMLETLC VQRLLSNMFF MHQLPPEHPD MQDLIQQNPE VSHLLDNSEI LCQTLELARH LAIQEIMQI QQPAQNPEHT LNPQPYLGLT TIPNGNIVLG QGYDNFNDHM LNGVHVHDDL EDDSFTALLA EQVLEQVQTP SLCQPPPQEQ WGQLSSSQVI FANSCGFLT TTEANLNTTN NASRENPAIA ATQGQSNVCA VQQPTEIPSL CTIPVIQEP EGDKDTTIPQ GSSAQWLEED LQQSEEQSSS QITRSMIQLL RTHPHMAAEM LLLLNPPQLN EQWRQQTATP LQSSQLSDLL LALANPKTSQ ALLQIEHGLQ LLATEAPALL PYIEPYLWGL GWLIPSSCDY LDTVPWTWNV QDTVEPQCPE SCHKSGTVLQ REQPLSGDSS HSLQAPEVRF SKEMECLQAM GFVNYNANLQ ALIATDGDTN AAIHKLKASQ GF
Source	Mammalian cell
Target Names	Ubqlnl
Protein Names	Recommended name: Ubiquilin-like protein
Expression Region	1-612
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.