



# Recombinant Human F-box/WD repeat-containing protein 9 (FBXW9)

<b>Product Code</b>	CSB-BP716858HU
<b>Abbreviation</b>	FBXW9
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5XUX1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MELPLGRCD D SRTWDDDDSDP ESETDPDAQA KAYVARVLSP PKSGLAFSRP SQLSTPAASP SASEPRAASR VSAVSEPGLL SLPELLEI CSYLDARLVL HVLSRVCHAL RDLVSDHVTW RLRALRRVRA PYPVVEEKNF DWPAACIALE QHLSRWAEDG RWVEYFCLAE GHVASVDSVL LLQGGSLCLS GSRDRNVNLW DLRQLGTESN QVLIKTLGTK RNSTHEGWVW SLAAQDHRVC SGSWDSTVKL WDMAADGQQF GEIKASSAVL CLSYLPDILV TGTYDKKVTI YDPRDRMETR DALMGWGGPS RGS DIPPHRP ITHEAGPALL KHQQLHSRPV LTL LADDRHI ISGSEDHTLV VVDRRANSVL QRLQLDSYLL CMSYQEPQLW AGDNQGLLHV FANRNGCFQL IRSFDVGHFSF PITGIQYSVG ALYTTSTDKT IRVHVPTDPP RTICTRRHDN GLNRVCAEGN LVVAGSGDLS LEVWRLQA
<b>Source</b>	Baculovirus
<b>Target Names</b>	FBXW9
<b>Protein Names</b>	Recommended name: F-box/WD repeat-containing protein 9 Alternative name(s): F-box and WD-40 domain-containing protein 9
<b>Expression Region</b>	1-488
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.