



Recombinant Human F-box/WD repeat-containing protein 9 (FBXW9)

Product Code	CSB-EP716858HU
Abbreviation	FBXW9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XUX1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MELPLGRCD D SRTWDDDDSDP ESETDPDAQA KAYVARVLSP PKSGLAFSRP SQLSTPAASP SASEPRAASR VSAVSEPGLL SLPELLEI CSYLDARLVL HVLSRVCHAL RDLVSDHVTW RLRALRRVRA PYPVVEEKNF DWPAACIALE QHLSRWAEDG RWVEYFCLAE GHVASVDSVL LLQGGSLCLS GSRDRNVNLW DLRQLGTESN QVLIKTLGTK RNSTHEGWVW SLAAQDHRVC SGSWDSTVKL WDMAADGQQF GEIKASSAVL CLSYLPDILV TGTYDKKVTI YDPRDRMETR DALMGWGGPS RGS DIPPHRP ITHEAGPALL KHQQLHSRPV LTL LADDRHI ISGSEDHTLV VVDRRANSVL QRLQLDSYLL CMSYQEPQLW AGDNQGLLHV FANRNGCFQL IRSFDVGHFSF PITGIQYSVG ALYTTSTDKT IRVHVPTDPP RTICTRRHDN GLNRVCAEGN LVVAGSGDLS LEVWRLQA
Source	E.coli
Target Names	FBXW9
Protein Names	Recommended name: F-box/WD repeat-containing protein 9 Alternative name(s): F-box and WD-40 domain-containing protein 9
Expression Region	1-488
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.