



Recombinant Nocardia farcinica Chaperone protein DnaJ 2 (dnaJ2)

Product Code	CSB-MP716879NAAA
Abbreviation	dnaJ2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5YNI2
Product Type	Recombinant Protein
Immunogen Species	Nocardia farcinica (strain IFM 10152)
Purity	≥85% (SDS-PAGE)
Sequence	MSQREWIEKD YYKELGVSST ASQDEIKKAY RKLARDLHPD ANPGDSKAE RFAVSEHA VLSDPAKRKE YDETRRLFAG GGYPRGGFGP GGAGGFSQGF DIGDIFGGTA AGDGGLGDIL GGLFNRRGTR TTSRPRRGAD VETETTLGFR EAAQGVTVPL RMTSPSPCTT CHGSGAKPGT SPRVCPICNG TGVVSRNQGA FGFSEPCDGC RGTGSIIDDP CVDCHGSGIQ NRTRTITVRI PPGVSDGQRI RLAGQGEAGL RGAPSGDLYV TVHVSQDKVF GRNGDDTLV LPVSYAELVL GTTVSVPTLD GRVGKVPVPPG TADGRTLVRV GRGVPKRGGG AGDLLVTVKV AVPQKLDDDA LEALERYREA EKASGFDPRA GWAGA
Source	Mammalian cell
Target Names	dnaJ2
Protein Names	Recommended name: Chaperone protein DnaJ 2
Expression Region	1-385
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.