



Recombinant Nocardia farcinica Chorismate synthase (aroC)

Product Code	CSB-EP716899NAAA
Abbreviation	aroC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5YTK8
Product Type	Recombinant Protein
Immunogen Species	Nocardia farcinica (strain IFM 10152)
Purity	≥85% (SDS-PAGE)
Sequence	MEAVLRWITA GESHGPALVA ILEGMVAGVE VTSEDISAQL ARRRRLGYGRG ARMKFEADKV TMVGGVRHGR TMGGPVAIEV ANSEWPKWTT VMSADPVDPA ELADLARNAP LTRPRPGHAD YSGMLKYNFD DARNVLERAS ARETAARVAA GTVARNFLRQ ALGVEVVS HV VAIGTAANTT GVVPTAADLA TVDESPVRAF DAAAEAAAMIA EIEAAKKDGD TLGGVVEVIV EGLPVGLGSF TSGENRLDSR LAAALMGIQA IKGVEVG DGF ETARRRGSQA HDEM R PGP DG VLRSTNRAGG LEGGMTNGEP LRVRAAMKPI STVPRALATV DMTTGEEAVA IHQRSDVCAV PAAGVVAEAM VALVVAQAAL EKFGGDSLTE TVDNITSYVK RISTRPHVPA DSRPQ
Source	E.coli
Target Names	aroC
Protein Names	Recommended name: Chorismate synthase EC= 4.2.3.5 Alternative name(s): 5-enolpyruvylshikimate-3-phosphate phospholyase
Expression Region	1-405
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.