



Recombinant Chicken Corepressor interacting with RBPJ 1 (CIR1)

Product Code	CSB-MP716936CH
Abbreviation	CIR1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZI03
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MGKSFANFMC KKDFHPASKS NIKKVWMAEQ KISYDKKKQE ELMQQYLKEQ ESYDNRLLMG DERVKGNLNF MYEAPPGAKK EEAKEQEGET EYKFEWQKVA PREKYAKDDM NIRDQPFGIQ VRNVRCIKCH KWGHVNTDRE CPLFGLSGIN ASSVSSDGSG PSMHPSELIA EMRNSGFALK QNVLGRNLTA NDPSQEFVAS EGDDDPEVEF LKSLSTKQKQ KLLRKLDRLE KKKKKKDRKK KKQQKKKGKS KHKKHKMSSS SSSSSETSGS SSDSETDSRD KAAQKKMYSK KRKKGKFSEA SSSDSEGKAK TRKEKLYEDL SSSHGNKDKD REKYRLLKQD SSAENSKRSS CEGERKSKSR YQRETERKER DSRQSADER RGRSPCTSHS SRRQRQTRKS PSQSPGEERH RKNETRSPST DGHRKEKKYR EGNGDRCSNV EYRERSRRSR SRSRERKKKR
Source	Mammalian cell
Target Names	CIR1
Protein Names	Recommended name: Corepressor interacting with RBPJ 1 Alternative name(s): CBF1-interacting corepressor
Expression Region	1-460
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.