



# Recombinant Chicken Snurportin-1 (SNUPN)

<b>Product Code</b>	CSB-EP716938CH-B
<b>Abbreviation</b>	SNUPN
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5ZI43
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Gallus gallus (Chicken)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MEELCAALAG GVALAAPNSP AAPHRLSAY KGRGDRLGQA ERRRLLCLQ RERRLDYVNH ARRLAEGDWA GVESEDEGGG DGDGEEEMEV DAGRRLPKRY ANQLMLSEWL VDVPVDLEQE WIVVVCVPGK RALVVASRGS TAAYTKSGFC VNRFPSELLPG GNRHNTMNEK VYCILDCIYN EAEQTYIILD VMCWRGHPVY DCQTDFFRFW LSSKIQEEEG LGEKSRINPY KFGVGLQNFPC TSESLCEVLT TNFPFEVDGL LFYHKQTHYT PGSTPLVGWL RPYMVPDILG LTVPATPLTA KPDYAGRQLQ QIIESKRSKK LAAGKAQPSA EAAARNGHYE LEHLSTPQPA NSAQQQEEAG SQMEN
<b>Source</b>	E.coli
<b>Target Names</b>	SNUPN
<b>Protein Names</b>	Recommended name: Snurportin-1 Alternative name(s): RNA U transporter 1
<b>Expression Region</b>	1-365
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	The nuclear import of the spliceosomal snRNPs U1, U2, U4 and U5, is dependent on the presence of a complex nuclear localization signal. The latter is composed of the 5'-2,2,7-terminal trimethylguanosine (m3G) cap structure of the U snRNA and the Sm core domain. This protein interacts specifically with m3G-cap and functions as an snRNP-specific nuclear import receptor. Alternatively spliced transcript variants encoding the same protein have been identified for this gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.