



# Recombinant Chicken CCR4-NOT transcription complex subunit 7 (CNOT7)

<b>Product Code</b>	CSB-EP716961CH-B
<b>Abbreviation</b>	CNOT7
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5ZJV9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Gallus gallus (Chicken)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MPAATVDHSQ RICEVWACNL DEEMKKIRQV IRKYNYVAMD TEFPGVVARP IGEFRSNADY QYQLLR CNVD LLKIIQLGLT FMNEQGEYPP GTSTWQFNFK FNLTEDMYAQ DSIELLTSG IQFKKHEEEG IETQYFAELL MTSGVVLCEG VKWLSFHSGY DFGYLIKILT NSNLPEEELD FFEILRLFFP VIYDVKYLK SCKNLKGGGLQ EVAEQLELER IGPQHQAGSD SLLTGMAFFK MREMFFEDHI DDAKYCGHLY GLGSGSSYVQ NGTGNAYEEE ANKQS
<b>Source</b>	E.coli
<b>Target Names</b>	CNOT7
<b>Protein Names</b>	Recommended name: CCR4-NOT transcription complex subunit 7 Alternative name(s): CCR4-associated factor 1 Short name= CAF-1
<b>Expression Region</b>	1-285
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This protein binds to an anti-proliferative protein, B-cell translocation protein 1, which negatively regulates cell proliferation. Binding of the two proteins, which is driven by phosphorylation of the anti-proliferative protein, causes signaling events in cell division that lead to changes in cell proliferation associated with cell-cell contact. The protein has both mouse and yeast orthologs. Alternate splicing of this gene results in two transcript variants encoding different isoforms.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.