



Recombinant Chicken Diphthamide biosynthesis protein 2 (DPH2)

Product Code	CSB-MP716971CH
Abbreviation	DPH2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZKI2
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MAAAFSSDGE AVLRRRLDPA AAAPRGDKDE FYEVDRAAAF VRDGGFRKVA LQFPDALLAD AAAVAARMEE VTGAEMYVLG DTTYGSCCVD EVAAEHVSAG AVVHYGPACL SPCRKLPLVH VFGRQPLDVG RCAEVFRELY PERQSRVVVL SDVVYAHAMG ELEKQLCHEY PNIIFSEVVC GDAPSPTLPG EVRQFGRRFH MEAAEELQDC SMFYVGAEGL ALTSFMLTNW RFPFSSFDPA TGHGRRETLN VNRALMRRLY LVERARDAHV VGILVGTGTV AGYLDVLEHL HQLVRRAGKR SYTLSVGKPN PAKLANFLEV DIFVLVACAQ NSLLDSSEFY RPIVTPYELE LACNPAREWT GNYLTDFRDL LPGACAHIEL PPAVPAAEAI PDVSLITGEM RATHLCDPLA PQPPSSTTLA CRDQTRALAE MSPAATFLES RSWRGLEQQL GKTAVSKAVQ GRRGIAIAYE DEGREQS
Source	Mammalian cell
Target Names	DPH2
Protein Names	Recommended name: Diphthamide biosynthesis protein 2 Alternative name(s): DPH2 homolog
Expression Region	1-477
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene is one of two human genes similar to the yeast gene dph2. The yeast gene was identified by its ability to complement a diphthamide mutant strain, and thus probably functions in diphthamide biosynthesis. Diphthamide is a post-translationally modified histidine residue present in elongation factor 2 (EF2) that is the target of diphtheria toxin ADP-ribosylation. Two transcript variants encoding different isoforms have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.