



# Recombinant Chicken Hsc70-interacting protein (ST13)

<b>Product Code</b>	CSB-BP716982CH
<b>Abbreviation</b>	ST13
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5ZLF0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Gallus gallus (Chicken)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MDSRKLGELR AFVRLCKQNP GLLHTEELGF LREWVESMGG TIPPAPASTS TDETSKGKAE EQPEEPVKSP EPESEESDLE IDNEGVIEPD NDDPQEMGDE NVEVTEEMMD QANEKKMEAI NALSEGDLQK AVNLFTDAIK LNPCLAILYA KRASVFKLQ KPNAAIRDCD RAIKINPDSA QTYKWRGKAH RLLGHWEEAA HDLALACKLD YEDASAMLK EVQPRAQKIA EHRRKYERKR EEKEIKERME RVKKAREEHE RAQREEEARR QAGGAQFGGF PGGFPGGFPG AMPGGMPGMA GMPGLNEILS DPEVLAAMQD PEVMAAFQDV AQNAPANMSKY QNNPKVMSLI TKLSAKFGSK P
<b>Source</b>	Baculovirus
<b>Target Names</b>	ST13
<b>Protein Names</b>	Recommended name: Hsc70-interacting protein Short name= Hip Alternative name(s): Protein FAM10A1 Protein ST13 homolog
<b>Expression Region</b>	1-361
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This protein is an adaptor protein that mediates the association of the heat shock proteins HSP70 and HSP90. This protein has been shown to be involved in the assembly process of glucocorticoid receptor, which requires the assistance of multiple molecular chaperones. The expression of this gene is reported to be downregulated in colorectal carcinoma tissue suggesting that is a candidate tumor suppressor gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.