



Recombinant Chicken Probable methylthioribulose-1-phosphate dehydratase (APIP)

Product Code	CSB-BP716988CH
Abbreviation	APIP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZLP2
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MAAASGHGLE LANGGDATQD KLHPRNLIPE LCRLFYGLGW VTGTGGGISL KHGNEIYIAP SGVQKERIQP EDMFVCDMNE QDISGPPAHK KLKKSQCTPL FMNAYTMRGA GAVIHTHska AVMATLLYPG NEFTITHQEM IKGIQKCSSG GYRYDDTLV VPIIENTPEE KDLKERMAKA MEEYPDSCAV LVRRHGVYVW GETWEKAKTM CECYDYLFDI AVQMKQHGLD PSKHPAGENG IL
Source	Baculovirus
Target Names	APIP
Protein Names	Recommended name: Probable methylthioribulose-1-phosphate dehydratase Short name= MTRu-1-P dehydratase EC= 4.2.1.109 Alternative name(s): APAF1-interacting protein homolog
Expression Region	1-242
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.