



Recombinant Chicken 14-3-3 protein beta/alpha (YWHAB)

Product Code	CSB-MP716989CH
Abbreviation	YWHAB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZLQ6
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MDKSELVQKA KLAEQAERYD DMAAAMKAVT EQGHELSENEE RNLLSVAYKN VVGARRSSWR VISSIEQKTE RNEKKQQMGR EYREKIEAEL QDICNDVLEL LDKYLIVNAT QPESKVFYLYL MKGDYYRYLS EVASGDNKQT TVANSQQAYQ EAFEISKKEM QPTHPIRLGL ALNFSVFYIE ILNSPEKACN LAKTAFDEAI AELDTLNEES YKDSTLIMQL LRDNLTLWTS ENQGDEGDAG EGEN
Source	Mammalian cell
Target Names	YWHAB
Protein Names	Recommended name: 14-3-3 protein beta/alpha
Expression Region	1-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a protein belonging to the 14-3-3 family of proteins, members of which mediate signal transduction by binding to phosphoserine-containing proteins. This highly conserved protein family is found in both plants and mammals. The encoded protein has been shown to interact with RAF1 and CDC25 phosphatases, suggesting that it may play a role in linking mitogenic signaling and the cell cycle machinery. Two transcript variants, which encode the same protein, have been identified for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.