



Recombinant Methanosarcina mazei V-type ATP synthase subunit D (atpD)

Product Code	CSB-MP717055MSN
Abbreviation	atpD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60188
Product Type	Recombinant Protein
Immunogen Species	Methanosarcina mazei (strain ATCC BAA-159 / DSM 3647 / Goe1 / Go1 / JCM 11833 / OCM 88) (Methanosarcina frisia)
Purity	≥85% (SDS-PAGE)
Sequence	MAQQDVKPTR SELINLKKKI KLSSESGHKLL KMKRDGLILE FFKILNEARN VRTELDAAFA KSTEKINLAS AVNGMVAVRS TAFTAKESPE IQLSGHNIMG VVVPKISSTG VRKSLYERGY GIIGTNSYID ETADAYEDLV EKIITAAELE TTMKRLLDEI EKTkRRVNAL EFKVPELID TMKYIRFMLE EMERENTFRL KRVKARMRD
Source	Mammalian cell
Target Names	atpD
Protein Names	Recommended name: V-type ATP synthase subunit D Alternative name(s): V-ATPase subunit D
Expression Region	1-209
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.