



Recombinant Mycoplasma hyopneumoniae tRNA pseudouridine synthase B (truB)

Product Code	CSB-BP717060MWK
Abbreviation	truB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q601U7
Product Type	Recombinant Protein
Immunogen Species	Mycoplasma hyopneumoniae (strain 232)
Purity	≥85% (SDS-PAGE)
Sequence	MADTLDLGSS AFGCKGSSPF LRTTMIGQPC VDFFFFFFF FKKKKKKMIT FLYKPKKISS AKFLRQWSKT NLIKKAGHAG TLDPLASGLL LVATEDDTKL LQYLDQKTKT YLAKIQFGFW STTYDAEGQI YGVEQPIKVT KDNLEQALNR LSESQKQVPP VFSSKKVSGK SAYHYARQGK QIELKPISIK ISKTILINF EKLQNCVIMW QVSRGCYIRS LADDLGKMLK TRAYLSDLER TKIGNFDKKF LNQSLKPQDL FDLQQVKLDL ENLELLQGK KINYFAKNSE LNTLIFKDEV VGFGKIINNV LITKKLFGNR IKKIINT
Source	Baculovirus
Target Names	truB
Protein Names	Recommended name: tRNA pseudouridine synthase B EC= 5.4.99.25 Alternative name(s): tRNA pseudouridine(55) synthase Short name= Psi55 synthase tRNA pseudouridylate synthase tRNA-uridine isomerase
Expression Region	1-327
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.