



Recombinant Aspartate carbamoyltransferase (pyrB)

Product Code	CSB-MP717063LAX
Abbreviation	pyrB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60252
Product Type	Recombinant Protein
Immunogen Species	Lactobacillus leichmannii
Purity	≥85% (SDS-PAGE)
Sequence	MTTVSQNQVV VETEDKQDNL LRLPYFVSVE QLSADDVLHL LQRAQYFKNG GEVPALSRPI FCTNMFFENS TRHTSFEVA ERLGLTVIP FDPSSHSSVNK GENLYDELT MASLGIELSV IRHPENAYYN EIIRPKEGQH LQMGLVNAGD GSGQHPSQSM LDMMTIYNEF GHFDGLKIMI VGDLTNSRVA RSNMEILNTL GAEVYFSGPE YWYNAEEFSK YGTYVKNIDD EIPELDVLML LRVQHERHNG AEAKTEQLFD AKDYNAAYGL NQRRYDMLKD DAIIMHPGPI NRGVEWDGDL VEAPKSRYAV QMHNGVFVRM AMIEAVLRGR KLGGLLE
Source	Mammalian cell
Target Names	pyrB
Protein Names	Recommended name: Aspartate carbamoyltransferase EC= 2.1.3.2 Alternative name(s): Aspartate transcarbamylase Short name= ATCase
Expression Region	1-336
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.