



Recombinant Methylococcus capsulatus

Glycerol-3-phosphate dehydrogenase [NAD (P)+] (gpsA)

Product Code	CSB-EP717098MFM-B
Abbreviation	gpsA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q604K1
Product Type	Recombinant Protein
Immunogen Species	Methylococcus capsulatus (strain ATCC 33009 / NCIMB 11132 / Bath)
Purity	≥85% (SDS-PAGE)
Sequence	MRHCISVLGA GSWG TALALL IARNGVPTLL WSHRREHAEA MAAERRNAHY LPHATFPDTL SVTPDLAEAA ATASDAVLVA VPSHAFRATL DLCAPHLQPE AGLMWATKGL DRTRGCLLHE TAREILGPDR ALAVLSGPTF AAEVVAELPS AITVGASSEA FGAKIVSWLH NAYFRAYTTD DLIGAQLGGA CKNVLAIAAG ISDGLGFGAN ARAALITRGL AEMMRLSLAL GARTETLMGL AGVGDLVLTA TDNQSRNRRF GLGVGKGRPR EEIHAEIGQE IEGILAAarli HELAGRKGVD MPITAQTYRV LYEGLPPEEA VRNLLLRDPK PEWDR
Source	E.coli
Target Names	gpsA
Protein Names	Recommended name: Glycerol-3-phosphate dehydrogenase [NAD(P)+] EC=1.1.1.94 Alternative name(s): NAD(P)H-dependent glycerol-3-phosphate dehydrogenase
Expression Region	1-335
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.