



Recombinant Mouse Melanocyte protein PMEL (Pmel), partial, Yeast

Product Code	CSB-YP717120MO
Abbreviation	Pmel
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60696
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Source	Yeast
Target Names	Pmel
Protein Names	Recommended name: Melanocyte protein PMEL Alternative name(s): Melanocyte protein Pmel 17 Premelanosome protein Silver locus protein Cleaved into the following 2 chains: 1. M-alpha 2. M-beta
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal His-tagged Tag-Free
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.