



Recombinant Mouse RAC-beta serine/threonine-protein kinase (Akt2)

Product Code	CSB-MP717138MO
Abbreviation	Akt2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60823
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MNEVSVIKEG WLHKRGEYIK TWRPRYFLLK SDGSFIGYKE RPEAPDQTLPLN PLNFSVAEC QLMKTERPRP NTFVIRCLQW TTVIERTFHV DSPDEREEWMM RAIQMVANSL KQRGPGEDAM DYKCGSPSDS STSEMMEVAV NKARAKVTMN DFDYLKLLGK GTFGKVLVR EKATGRYYAM KILRKEVIA KDEVAHTVTE SRVLQNTRHP FLTALKYAFQ THDRLCFVME YANGGELFFH LSRERVFTED RARFYGAEIV SALEYLHSRD VVYRDIKLEN LMLDKDGHK ITDFGLCKEG ISDGATMKTF CGTPEYLAPE VLEDNDYGRA VDWWGLGVVM YEMMCGRLPF YNQDHERLFE LILMEEIRFP RTLGPPEAKSL LAGLLKDPK QRLGGGPSDA KEVMEHRFFL SINWQDVVQK KLLPPFKPQV TSEVDTRYFD DEFTAQSITI TPPDRYDSL DPLELDQRTHF PQFSYSASIR E
Source	Mammalian cell
Target Names	Akt2
Protein Names	Recommended name: RAC-beta serine/threonine-protein kinase EC= 2.7.11.1 Alternative name(s): Protein kinase Akt-2 Protein kinase B beta Short name= PKB beta RAC-PK-beta
Expression Region	1-481
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene is a putative oncogene encoding a protein belonging to a subfamily of serine/threonine kinases containing SH2-like (Src homology 2-like) domains. The gene was shown to be amplified and overexpressed in 2 of 8 ovarian carcinoma cell lines and 2 of 15 primary ovarian tumors. Overexpression contributes to the malignant phenotype of a subset of human ductal pancreatic cancers. The encoded protein is a general protein kinase capable of phosphorylating several known proteins.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.