



Recombinant Mouse ELAV-like protein 2 (Elavl2)

Product Code	CSB-MP717142MO
Abbreviation	Elavl2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60899
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	METQLSNGPT CNNTANGPTT VNNNCSSPVD SGNTEDSKTN LIVNYLPQNM TQEELKSLFG SIGEIESCKL VRDKITGQSL GYGfVNYIDP KDAEKAINTL NGLRLQTKTI KVSYPARPSSA SIRDANLYVS GLPKTMTQKE LEQLFSQYGR IITSRLVDQ VTGISRGVGF IRFDKRIEAE EAIKGLNGQK PPGATEPITV KFANNPSQKT NQAILSQLYQ SPNRRYPGPL AQQAQRFRLD NLLNMAYGVK SRFSPMTIDG MTSLAGINIP GHPGTGWCIF VYNLAPDADE SILWQMFGPF GAVTNVKVIR DFNTNKCKGF GFVTMTNYDE AAMAIASLNG YRLGDRVLQV SFKTNKTHKA
Source	Mammalian cell
Target Names	Elavl2
Protein Names	Recommended name: ELAV-like protein 2 Alternative name(s): ELAV-like neuronal protein 1 Hu-antigen B Short name= HuB Nervous system-specific RNA-binding protein Mel-N1
Expression Region	1-360
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.