



# Recombinant Mouse Protein p13 MTCP-1 (Mtcp1)

<b>Product Code</b>	CSB-BP717148MO
<b>Abbreviation</b>	Mtcp1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q60945
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAREDVGAPP DHLWVHQEGV YRDEYQRTWV AVVEEETSFL KARVQQVQVP LGDATKPSHL LTSQLPLMWQ LYPEERYMDN NSRLWQIQHH LMVRGVQELL LKLLPDD
<b>Source</b>	Baculovirus
<b>Target Names</b>	Mtcp1
<b>Protein Names</b>	Recommended name: Protein p13 MTCP-1 Short name= p13MTCP1 Alternative name(s): Mature T-cell proliferation-1 type B1 Short name= MTCP-1 type B1
<b>Expression Region</b>	1-107
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene was identified by involvement in some t(X;14) translocations associated with mature T-cell proliferations. This region has a complex gene structure, with a common promoter and 5 exon spliced to two different sets of 3 exons that encode two different proteins. This gene represents the upstream 13 kDa protein that is a member of the TCL1 family. This protein may be involved in leukemogenesis.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.