



# Recombinant *Methylococcus capsulatus* Chaperone surA (surA)

<b>Product Code</b>	CSB-EP717161MFM-B
<b>Abbreviation</b>	surA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q60B78
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Methylococcus capsulatus</i> (strain ATCC 33009 / NCIMB 11132 / Bath)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	AV DRIVAVVDDG VILESELVRK VDEIKRSLRA SRASLPPDSV LVRQVLERMI VDKIQIQMAE KMGIQVDDDT LRMAVSQIAQ RNNLTPDQFR RSLAREGIDY GDFLDQVRSE IAMGRLRASQ INNQIKISDR EVEHYLEAQQ GSGAVADREY RLGHILIATP REASPDEVKK ARERADR VVK ELKAGLDFKD ASIRYSDDPQ ALEGGDLGWR KLSEIPSHIA EVVGGMKDGE VSDPIRSPGG YHIVKMLAMR GVGEAKLTKT HVRHILIRPN EVLSDEDAKN KLLALKTRIE NGDDFAELAR GHSDDKGS AI KGGDLGWVKP GALVPPFEEA MNALDENQLS DPVQTQFGWH LIQVLERQES SDTNEVLKNR ARDELFKRKV DEETELWLRK IRDEAYVEIR LDETPASPG E DAPAGEDSPE TFMR
<b>Source</b>	E.coli
<b>Target Names</b>	surA
<b>Protein Names</b>	Recommended name: Chaperone surA Alternative name(s): Peptidyl-prolyl cis-trans isomerase surA Short name= PPlase surA EC= 5.2.1.8 Rotamase surA
<b>Expression Region</b>	29-454
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.