



Recombinant Methylococcus capsulatus

Chaperone surA (surA)

Product Code	CSB-YP717161MFM
Abbreviation	surA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60B78
Product Type	Recombinant Protein
Immunogen Species	Methylococcus capsulatus (strain ATCC 33009 / NCIMB 11132 / Bath)
Purity	≥85% (SDS-PAGE)
Sequence	AV DRIVAVVDDG VILESELVRK VDEIKRSLRA SRASLPPDSV LVRQVLERMI VDKIQIQMAE KMGIVDDDDT LRMAVSQIAQ RNNLTPDQFR RSLAREGIDY GDFLDQVRSE IAMGRLRASQ INNQIKISDR EVEHYLEAQG GSGAVADREY RLGHILIATP REASPDEVKK ARERADRVVK ELKAGLDFKD ASIRYSDDPQ ALEGGDLGWR KLSEIPSHIA EVVGGMKDGE VSDPIRSPGG YHIVKMLAMR GVGEAKLTKT HVRHILIRPN EVLSDEDAKN KLLALKTRIE NGDDFAELAR GHSDDKGSAL KGGDLGWVKP GALVPPFEEA MNALDENQLS DPVQTQFGWH LIQVLERQES SDTNEVLKNR ARDELFKRKV DEETELWLRK IRDEAYVEIR LDETPASPGE DAPAGEDSPE TFMR
Source	Yeast
Target Names	surA
Protein Names	Recommended name: Chaperone surA Alternative name(s): Peptidyl-prolyl cis-trans isomerase surA Short name= PPlase surA EC= 5.2.1.8 Rotamase surA
Expression Region	29-454
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.