



Recombinant Methylococcus capsulatus S-adenosylmethionine synthase (metK)

Product Code	CSB-BP717166MFM
Abbreviation	metK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60CG7
Product Type	Recombinant Protein
Immunogen Species	Methylococcus capsulatus (strain ATCC 33009 / NCIMB 11132 / Bath)
Purity	≥85% (SDS-PAGE)
Sequence	MSNSFIFTSE SVSEGHPDKI ADQVSDAVLD AILAQDPKAR VAVETLVKTG MVLAGEVTT SAWVDTEELV RKVVHEIGYD NPEIGFDWQS CAVLTAIGKQ SPDIAMGVDE FDDHEQGAGD QGLMFGYATN ETDVLMPIPI TYAHLVQRQ AEVRKHKVLP WLRPDAKSQV TFRYEEGKPV AVDAVVLSTQ HAPDVAQKDI REAVMEEIY HVLPREWLHK DTKYFINPTG SFVIGGPVGD CGLTGRKIIV DTYGGMARHG GGAFSGKDPS KVDRSAAYMG RYVAKNIVAA GLAERCEIQV SYAIGVAEPT SISIETFGTG RIAEETLVKL VREHFDLRPK GLITMLDLLR PIYRPTAAYG HFGRTGETFT WERTDKADLL REAAGLPPIA
Source	Baculovirus
Target Names	metK
Protein Names	Recommended name: S-adenosylmethionine synthase Short name= AdoMet synthase EC= 2.5.1.6 Alternative name(s): MAT Methionine adenosyltransferase
Expression Region	1-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.