



Recombinant *Oryza sativa* subsp. japonica Protein N-terminal glutamine amidohydrolase (Os05g0387900, LOC_Os05g32190)

Product Code	CSB-MP717172OFG
Abbreviation	Os05g0387900, LOC_Os05g32190
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60E61
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MADDRVAGGA T P P P P P P P P P L D A S A F T H T P Y Y C E E N V H L L C K E L I R S G I S D P A G T N L Y A V F I S N E E K K V P L W Y Q K A S H S G D G F V L W D Y H V I C I Q S R R K N G E V L D L V W D L D S S L P F P C S F I Q Y V S D A I R P L S F G N S T Y R R L F R V I H A P V F L R S F A S D R S H M K D H A G N W I Q L P P K Y E S I V A E D G T T N N L N E Y I T M S M D D V K D L E S M A D D V Y S S K H G V V I N E T I L P E F F S R L P G
Source	Mammalian cell
Target Names	Os05g0387900
Protein Names	Recommended name: Protein N-terminal glutamine amidohydrolase EC= 3.5.1.- Alternative name(s): Protein NH ₂ -terminal glutamine deamidase Short name= N-terminal Gln amidase Short name= Nt(Q)-amidase
Expression Region	1-231
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.