



Recombinant DOMON domain-containing protein CBG21755 (CBG21755)

Product Code	CSB-EP717184CXX
Abbreviation	CBG21755
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60QM6
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis briggsae
Purity	≥85% (SDS-PAGE)
Sequence	ATCDFENEVA SMSWMVKGDV LQINFEHKNL TENRWASIAF GNGPGMTNLE SIIFSRDDNN VITTNTGFTP KKKKVVVDDV SYVTVRNVQL KGDLLRITVT RPLGPAGPRD FSLDQCVNWI VPPGGALSNG KFKKHHGQIY FVKDVCAAKC TAQRMMRVLS NRIH
Source	E.coli
Target Names	CBG21755
Protein Names	Recommended name: DOMON domain-containing protein CBG21755
Expression Region	21-184
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.