





Recombinant Rat Interleukin-17A (II17a)

Product Code	CSB-YP717213RA
Abbreviation	II17a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q61453
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	AVL IPQSSVCPNA EANNFLQNVK VNLKVINSLS SKASSRRPSD YLNRSTSPWT LSRNEDPDRY PSVIWEAQCR HQRCVNAEGK LDHHMNSVLI QQEILVLKRE PEKCPFTFRV EKMLVGVGCT CVSSIVRHAS
Source	Yeast
Target Names	II17a
Protein Names	Recommended name: Interleukin-17A Short name= IL-17 Short name= IL-17A Alternative name(s): Cytotoxic T-lymphocyte-associated antigen 8 Short name= CTLA-8
Expression Region	18-150
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is a proinflammatory cytokine produced by activated T cells. This cytokine regulates the activities of NF-kappaB and mitogen-activated protein kinases. This cytokine can stimulate the expression of IL6 and cyclooxygenase-2 (PTGS2/COX-2), as well as enhance the production of nitric oxide (NO). High levels of this cytokine are associated with several chronic inflammatory diseases including rheumatoid arthritis, psoriasis and multiple sclerosis.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL.We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



CUSABIO TECHNOLOGY LLC



of lyophilized form is 12 months at -20°C/-80°C.