



# Recombinant Mouse Probable ATP-dependent RNA helicase DDX4 (Ddx4)

<b>Product Code</b>	CSB-EP717216MO
<b>Abbreviation</b>	Ddx4
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q61496
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGDEDWEAEI LKPHVSSYVP VFEKDKYSSG ANGDTFNRTS ASSDIGESSK KENTSTTGGF GRGKGFGRG FLNNKFEEGD SSGFWKESNN DCEDNQTRSR GFSKRGGCQD GNDSEASGPF RRGGRGSFRG CRGGFGLGRP NSESDQDQGT QRGGGLFGSR KPAASDSGNG DTYQSRSGSG RGGYKGLNEE VVTGSGKNSW KSETEGGESS DSQGPKVTYI PPPPEDEDS IFAHYQTGIN FDKYDTILVE VSGHDAPPAI LTFEEANLCQ TLNNNIAKAG YTKLTPVQKY SIPIVLAGRD LMACAQTGSG KTA AFLLPIL AHMMRDGITA SRFKELQEPE CIIVAPTREL INQIYLEARK FSFGTCVRAV VIYGGTQFGH SVRQIVQGCN ILCATPGRLM DIIGKEKIGL KQVKYLVLDE ADRMLDMGFG PEMKKLISCP GMPKSKEQRQT LLFSATFPPEE IQRLAGDFLK SSYLFVAVGQ VGGACRDVQQ TILQVGQYSK REKLVEILRN IGDERTMVFV ETKKKADFA TFLCQEKIST TSIHGDREQR EREQALGDFR CGKCPVLVAT SVAARGLDIE NVQHVINFDL PSTIDEYVHR IGRTGRCGNT GRAISFFDTD SDNHLAQLPV KVLSDAQQDV PAWLEEIAFS TYVPPSFSSS TRGGAVFASV DTRKNYQGKH TLNTAGISSS QAPNPVDDDES WD
<b>Source</b>	E.coli
<b>Target Names</b>	Ddx4
<b>Protein Names</b>	Recommended name: Probable ATP-dependent RNA helicase DDX4 EC=3.6.4.13 Alternative name(s): DEAD box protein 4 Mvh Vasa homolog
<b>Expression Region</b>	1-702
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	DEAD box proteins, characterized by the conserved motif Asp-Glu-Ala-Asp (DEAD), are putative RNA helicases. They are implicated in a number of cellular processes involving alteration of RNA secondary structure such as translation



initiation, nuclear and mitochondrial splicing, and ribosome and spliceosome assembly. Based on their distribution patterns, some members of this family are believed to be involved in embryogenesis, spermatogenesis, and cellular growth and division. This gene encodes a DEAD box protein, which is a homolog of VASA proteins in *Drosophila* and several other species. The gene is specifically expressed in the germ cell lineage in both sexes and functions in germ cell development. Multiple transcript variants encoding different isoforms have been found for this gene.

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.