



Recombinant Mouse ELAV-like protein 4 (Elavl4)

Product Code	CSB-EP717224MO-B
Abbreviation	Elavl4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q61701
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MEWNLKMI STMEPQVSNP PTSNTSNGPS SNNRNCPSM QTGAATDDSK TNLIVNYLPQ NMTQEEFRSL FGSIGEIESC KLVRDKITGQ SLGYGFVNYI DPKDAEKAIN TLNGLRLQTK TIKVSYARPS SASIRDANLY VSGLPKTMTQ KELEQLFSQY GRIITSRILV DQVTGVSRGV GFIRFDKRIE AEEAIKGLNG QKPSGATEPI TVKFANNPSQ KSSQALLSQL YQSPNRRYPG PLHHQAQRFR LDNLLNMAYG VKRLMSGPVP PSACPPRFSP ITIDGMTSLV GMNIPGHTGT GWCIFVYNLS PDSDESVLWQ LFGPFGAVNN VKVIRDFNTN KCKGFGFVTM TNYDEAAMAI ASLNGYRLGD RVLQVSFKTN KAHKS
Source	E.coli
Target Names	Elavl4
Protein Names	Recommended name: ELAV-like protein 4 Alternative name(s): Hu-antigen D Short name= HuD Paraneoplastic encephalomyelitis antigen HuD
Expression Region	1-385
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.