



Recombinant Mouse Protein phosphatase 1 regulatory subunit 14B (Ppp1r14b)

Product Code	CSB-BP717249MO
Abbreviation	Ppp1r14b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q62084
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	ADSGPAGGA ALAAPAPGPG SGSTGPRVYF QSPPGAAGEG PGGADDDGPV RRQGKVTVKY DRKELRKRLN LEEWILEQLT RLYDCQEEI PELEIDVDEL LDMESDDTRA ARVKELLVDC YKPTEAFISG LLDKIRGMQK LSTPQKK
Source	Baculovirus
Target Names	Ppp1r14b
Protein Names	Recommended name: Protein phosphatase 1 regulatory subunit 14B Alternative name(s): Phosphatase holoenzyme inhibitor 1 Short name= PHI-1 Phospholipase C-beta-3 neighbouring gene protein
Expression Region	2-147
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.