



Recombinant Mouse Protein SSXT (Ss18)

Product Code	CSB-BP717257MO
Abbreviation	Ss18
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q62280
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	SVAFAAPRQ RGKGEITPAA IQKMLDENNH LIQCIMDYQN KGKASECSQY QQILHTNLVY LATIADSNQN MQSLLPAPPT QTMPMGPGGM SQSGPPPPPR SHNMPSDGMV GGGPPAPHMV NQMNGQMPGP NHMPMQGPGP SQLSMTNSSM NMPSSSHGSM GGYNHSPVSS QSMQVQNMNT MSQGQPMGNY GPRPNMNMQP NQGPMMHQQP PSQQYNMPPG GAQHYQGQQA PMGLMGQVNQ GSHMMGQRQM PPYRPPQQGP PQQYSGQEDY YGDQYSHGGQ GPPEGMNQQY YPDGHNDYGY QQPSYPEQGY DRPYEDSSQH YYEGGNSQYG QQQDAYQGPP PQQGYPPQQQ QYPGQQGYPG QQQSYGPSQG GPGPQYPNYP QGQGQQYGGY RPTQPGPPQP PQRPHYGYDQ GQYGNYYQ
Source	Baculovirus
Target Names	Ss18
Protein Names	Recommended name: Protein SSXT Alternative name(s): Protein SYT Synovial sarcoma-associated Ss18-alpha
Expression Region	2-418
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.