



Recombinant Mouse Endophilin-A3 (Sh3gl3)

Product Code	CSB-BP717263MO
Abbreviation	Sh3gl3
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q62421
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MSVAGLKKQF HKASQLFSEK ISGAEGTKLD EEFLNMEKKI DITSKAVAEI LSKATEYLQP NPAYRAKLG M LNTVSKLRGQ VKATGYPQTE GLLGDCMLKY GKELGEDSAF GNSLVDVGEA LKLMAEVKDS LDINVKQTFI DPLQLLQDKD LKEIGHHLRK LEGRRLDYDY KKRRVGIPIE EEIRQAVEKF EESKELAERS MFNFLENDVE QVSQLAVFVE AALDYHRQST EILQELQSKL ELRISLASKV PKREFMPKPV NMSSTDANGV GPSSSSKTPG TDTPADQPCC RGLYDFEPEN EGELGFKEGD IITLTNQIDE NWYEGMLRGE SGFFPINYVE VIVPLPP
Source	Baculovirus
Target Names	Sh3gl3
Protein Names	Recommended name: Endophilin-A3 Alternative name(s): Endophilin-3 SH3 domain protein 2C SH3 domain-containing GRB2-like protein 3 SH3p13
Expression Region	1-347
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>