



Recombinant Mouse Transducin-like enhancer protein 1 (Tle1)

Product Code	CSB-MP717264MO
Abbreviation	Tle1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q62440
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MFPQSRHPTPHQAAGQPFKFTIPESLDRIKEEFQFLQAQYHSLKLECEKLASE KTEMQRH YVMYYEMSYGLNIEMHKQTEIAKRLNTICAQVIPFLSQEHQQQVAQAVERAQK VTMAELN AIIGQQQLQAQHLSHGHGPPVPLTPHPSGLQPPGIPPLGGSASLLALSSALSG QSHLAIK DDKKHHDAERHRDRESGTSNSLLVPDSLSTDKRRNGPEFSSDIKKRKVDDK DNYDSDGD KSDDNLVVDVSNEDPSSPHASPTHSPRENGIDKNRLLKKDASGSPASTASSG SSSSLKSK EVSLHEKANTPVLKSSTPTPRSDMPTPGTSATPGLRPGLGKPPAMEPLVNQA AAGLRTP AVPGPYPAFPGMVPHAGMNGELTSPGAAYAGLHSMSPQMSAAAAAAAAAAV AYGRSPMVG FDPPPHMRVPSIPPNLAGIPGGKPAYSFHVTADGQMMPVFPFPPDALIGPGIPR HARQINT LNHGEVCAVTISNPTRHVYTGKGKCVKWDISHPGNKSPVSQDCLNRDNYI RSCKLLP DGCTLIVGGEASTLSIWDLAAPTPIKAELTSSAPACYALAIAPDSKVCFSCCSD GNIIV WDLHNQTLVRQFQGHTDGASCIDISNDGTLWTGGLDNTVRSWDLREGRQL QQHDFTSQI FSLGYCPTGEWLAVGMESSNVEVLHVKNPKDKYQLHLHESCVLSLKFAAYCGKW FVSTGKDN LLNAWRTPYASIFQSKESSSVLSCDISVDDKYIVTGSGDKKATVYEVY
Source	Mammalian cell
Target Names	Tle1
Protein Names	Recommended name: Transducin-like enhancer protein 1 Alternative name(s): Groucho-related protein 1 Short name= Grg-1



Expression Region	1-770
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.