



Recombinant Mouse Transducin-like enhancer protein 1 (Tle1)

Product Code	CSB-YP717264MO
Abbreviation	Tle1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q62440
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MFPQSRHPTPHQAAGQPFKFTIPESLDRIKEEFQFLQAQYHSLKLECEKLASE KTEMQRH YVMYYEMSYGLNIEMHKQTEIAKRLNTICAQVIPFLSQEHQQQVAQAVERAQKQ VTMAELN AIIGQQQLQAQHLSHGHGPPVPLTPHPSGLQPPGIPPLGGSASLLALSSALSG QSHLAIK DDKKHHDAERHRDRESGTSNSLLVPDSL RSTDKRRNGPEFSSDIKKRKVDDK DNYDSDGD KSDDNLVVDVSNEDPSSPHASPTHSPRENGIDKNRLLKKDASGSPASTASSG SSSSLKSK EVSLHEKANTPVLKSSSTPTPRSDMPTPGTSATPGLRPGLGKPPAMEPLVNQA AAGLRTP AVPGPYPAFPGMVPHAGMNGELTSPGAAYAGLHSMSPQMSAAAAAAAAAAV AYGRSPMVG FDPPPHMRVPSIPPNLAGIPGGKPAYSFHVTADGQM QPVFPFPDALIGPGIPR HARQINT LNHGEVCAVTISNPTRHVYTGKGKCVKWDISHPGNKSPVSQDCLNRDNYI RSCKLLP DGCTLIVGGEASTLSIWDLAAPT PRIKAELTSSAPACYALAI SPDSKVCFSCCSD GNI AV WDLHNQTLVRQFQGHTDGASCIDISNDG TKLWTGGLDNTVRSWDLREGRQL QQHDFTSQI FSLGYCPTGEWLAVGMESSNVEVLHV NKPDKYQLHLHESCVLSLK FAYCGKW FVSTGKDN LLNAWRTPY GASIFQSKESSSVLSCDISVDDKYIVTGSGDKKATVYEVY
Source	Yeast
Target Names	Tle1
Protein Names	Recommended name: Transducin-like enhancer protein 1 Alternative name(s): Groucho-related protein 1 Short name= Grg-1



Expression Region	1-770
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.