



Recombinant Mouse Cochlin (Coch)

Product Code	CSB-MP717267MO
Abbreviation	Coch
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q62507
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	EGAV PIPVTCFTRG LDIRKEKADV LCPGGCSLEE FSVFGNIVYA SVSSICGAAV HRGVIGTSGG PVRVYSLPGR ENYSSVDANG IQSQMLSRWS ASFAVTKGKS STQEATGRAV STAHPSPGKR LKKTPEKKTG NKDCKADIAF LIDGSFNIGQ RRFNLQKNFV GKVALMLGIG TEGPHVGLVQ ASEHPKIEFY LKNFTSAKDV LFAIKEVGFR GGNSTGKAL KHTAQKFFTA DTGVRKGIPK VVVVFIDGWP SDDIEEAGIV AREFGVNVFI VSVAKPIPEE LGMVQDVAFV DKAVCRNNGF FSYHMPNWFV TTKYVKPLVQ KLCTHEQMMC SKTCYNSVNI AFLIDGSSSV GDSNFRLMLE FVSNIAKTFE ISDIGAKIAA VQFTYDQRTE FSFTDYNTKE NVLAVLANIR YMSGGTATGD AIAFTVRNVF GPIRDSPNKN FLVIVTDGQS YDDVRGPAAA AHDAGITIFS VGVAWAPLDD LRDMSKPKE SHAFFTREFT GLEPIVSDVI RGICRDFLES QQ
Source	Mammalian cell
Target Names	Coch
Protein Names	Recommended name: Cochlin Alternative name(s): COCH-5B2
Expression Region	27-552
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is highly conserved in human, mouse, and chicken, showing 94% and 79% amino acid identity of human to mouse and chicken sequences, respectively. Hybridization to this gene was detected in spindle-shaped cells located along nerve fibers between the auditory ganglion and sensory epithelium. These cells accompany neurites at the habenula perforata, the opening through which neurites extend to innervate hair cells. This and the pattern of expression of this gene in chicken inner ear paralleled the histologic findings of acidophilic deposits, consistent with mucopolysaccharide ground substance, in temporal bones from DFNA9 (autosomal dominant nonsyndromic sensorineural deafness 9) patients. Mutations that cause DFNA9 have been



reported in this gene. Alternative splicing results in multiple transcript variants encoding the same protein. Additional splice variants encoding distinct isoforms have been described but their biological validities have not been demonstrated.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.