



Recombinant Rat NGFI-A-binding protein 1 (Nab1)

Product Code	CSB-MP717278RA
Abbreviation	Nab1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q62722
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MATALPRTLGL ELQLYRILQK ANLLSYFDAF IQQGGDDVQQ LCEAGEEEFL EIMALVGMAS KPLHVRRLQK ALRDWVTNPG LFNQPLTSLP VSSIPIYKLP EGSPTWLGIS CNSYERSSSA REPHLKVPKC AATTCVQSLG QGKSEVGSLA LQSVSESRLW QGHHTTESEH SLSPADLGSP ASPKESSEAL DAAAALSVAE CVERMAPTLP KSDLNEVKEL LKNNKKLAKM IGHIFEMSDE DPHKEEEIRK YSAIYGRFDS KRKDGKHLTL HELTVNEAAA QLCVKDNALL TRRDELFALA RQVSREVTYK YTYRTTRLKC GERDELSPKR IKMEDGFPDF QESVPTLFQQ ARAKSEELAG LSSQAEKGMA KQMELLCAQA GYERLQQERR LMAGLYRQSS GEQSPDGGLP SDGSDGQGER PLNLRIPSVQ NRQPHHFVVD GELSRLYSNE VKSHSSES LG ILKDYPHSAF TLEKKVIKTE PEDSR
Source	Mammalian cell
Target Names	Nab1
Protein Names	Recommended name: NGFI-A-binding protein 1 Alternative name(s): EGR-1-binding protein 1
Expression Region	1-485
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.